(90 × 210mm Size)

**FUJIFILM**

**Wako**

〈For Research Use Only〉 Code No. 296-63801 (1,300 Tests)

## LabAssay™ Phospholipid (Choline Oxidase·DAOS method)

For Research Use Only Not for Diagnostic Use

### [Introduction]

Phospholipids are a class of lipids, and a major component of all biological membranes, along with glycolipids, cholesterol and proteins.

LabAssay™ Phospholipid is based on an enzymatic methods using *N*-ethyl-*N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) as a blue pigment. This kit is used for the quantitative determination of phospholipids in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

### [Kit contents]

(1)	Buffer (50mmol/L Good' s buffer, pH7.5)	50mL	8vials
(2)	Chromogen Substrate (when reconstituted)	For 50mL	8vials
	Phospholipase D           0.47units/mL Choline Oxidase         2.0units/mL Peroxidase               4.2units/mL <i>N</i> -Ethyl- <i>N</i> -(2-hydroxy-3-sulfopropyl)-3,5- dimethoxyaniline sodium salt (DAOS) 0.77mmol/L 4-Aminoantipyrine       0.24mmol/L Ascorbate Oxidase       3.9units/mL		
(3)	Standard Solution (Choline Chloride           54mg/dL (Corresponding to 300mg/dL Phospholipid)	10mL	2vials

### [Materials and apparatuses to be prepared]

- 96-well microplate (transparent type)
- Micropipette
- Incubator maintaining at 37°C\*
- Plate mixer\*
- Microplate reader with a 600nm wavelength filter  
(\* if the microplate reader is not equipped.)

### [For test tube method]



- Test tube
- Pipette
- Spectrophotometer or colorimeter with 600 nm wavelength filter

### [Assay principle]


Phospholipids (lecithin, sphingomyelin and lysolecithin) in a sample are

—1/8—


hydrolyzed to choline in a reaction catalyzed by phospholipase D. Choline so formed is oxidized by choline oxidase in a reaction that produces hydrogen peroxide. The hydrogen peroxide produced causes DAOS and 4-Aminoantipyrine to undergo a quantitative oxidative condensation catalysed by peroxidase (POD), producing a blue pigment. The amount of phospholipids contained in the sample is determined by measuring the absorbance of the blue color.

Choline Oxidase

POD

### [Preparation of reagents to be used]

#### ① Color Reagent :

Prepare Color Reagent by Dissolving 1 vial of Chromogen substrate (for 50mL) to 50mL of Buffer. After reconstitution, the solution should be stored at 2-10°C and used within one week.

#### ② Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard.

No.	Phospholipid Standard	Distilled or deionized water	Sample volume	Phospholipid Concentration
1	100 μL	100 μL	2 μL	150 mg/dL
2	Undiluted	—	2 μL	300 mg/dL
3	Undiluted	—	4 μL	596.1 mg/dL* <sup>1</sup>

\* 1 The test sample volume is usually 2 μL, but 4 μL is taken in this case. The phospholipids concentration must be corrected accordingly as indicated in the table above.

### [Procedure]

#### (1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Sample	Serum 2 μL	Standard solution 2 μL	—
Color reagent	300 μL	300 μL	300 μL
	Mix well and incubate at 37°C for 5 min. Measure the absorbance at 600nm* <sup>2</sup> of the test sample and standard solution with the blank solution as the control.		

\*2 In two wavelength assay, measure using Main wavelength 600nm/Sub wavelength 700nm.

—2/8—


(2) Assay in a test tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Sample	Serum 0.02 mL	Standard solution 0.02 mL	—
Color reagent	3 mL	3 mL	3 mL
	Mix well and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600nm Spectrophotometer : 600nm*2		

\*2 In two wavelength assay, measure using Main wavelength 600nm/Sub wavelength 700nm.

[Standard curve] [assay in a microplate]

Microplate reader : SAFIRE (TECAN)

[Performance]

(1) Sensitivity [assay in a test tube]

- The absorbance is below 0.14, when measuring purified water as a sample.
- The absorbance is 0.12 to 0.58, when measuring 300mg/dL phospholipid as a sample.

(2) Specificity

- The phospholipid concentration is less than  $\pm 10\%$ , when measuring the known concentration of control serum as a sample.

[Usage Notes]

(1) Sample

- Ascorbic acid, bilirubin and hemolysis may not significantly affect the assay.
- This method measures lecithin, sphingomyelin, lysolecithin but not cephalin.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other

than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.

- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- The reaction completes in about 3 min. It is enough to incubate at 37 °C for 5 min. When the incubation is continued, the absorbance gradually decreases. Therefore, the incubation should be done within 15 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This kit is for research use only. Not for diagnostic use.

(3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette should be used with a safety pipette filler.

(4) Waste

- The waste should be processed appropriately according to the law.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

**Expiration date** : 12 months after the manufacture

**Storage** : Store at 2 ~ 10°C

**Package** : 1,300 tests

[Reference]

1. Takayama, M., Itoh, S., Nagasaki, T. and Tanimizu, I. : *Clin. Chim. Acta.*, **79**, 93 (1977).

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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